



ENAR ORDER FORM

19 MAY 21

Date _____ Referred by _____

Name or Co. _____

Address _____

Delivery _____

Work/Home Ph _____ Mob _____

Email _____



Enlightened Therapies Pty Ltd

ABN 96 099 251 117

info@enlightenedtherapies.com

30 Park St Coledale NSW 2515

PO Box 3130 Austinmer NSW 2515

Phone + 61 (2) 4268 2222

Toll Free 1300 305 107

www.enar.com.au



Therapy System



PAYMENT PLANS

ET = INTEREST FREE
REPAYMENT OPTIONS *

EG : Pay 50% upfront &
Pay Balance in 50 days

BY ARRANGEMENT *

FREE SUPPORT

- ✓ ENAR devices are sold with a Two Year Repair or Replacement Warranty.
- ✓ Training Resources incl DVD & Colour Manual
- ✓ Free local call cost
Help-Line 1300 305 107
- ✓ Free Life-Time Ongoing ENAR Advice & Support

OUR GUARANTEE

**Enlightened Therapies
Ongoing Support Guarantee**

ET will respond personally to you within 48 hours and will provide help on how best to use your ENAR therapy device and accessories.

ENAR two year replacement warranty is for manufacturing defects. (Not for user misuse.)

After its two year warranty period your device may be sent to us for assessment, service, repair, and / or replacement options.

YOUR FEEDBACK

What are you hoping to treat?

How did you hear about ET?

PRODUCTS	PRICE	QTY	TOTAL
ENAR Therapy System Device Special * Save \$100 *	\$890		
Special includes 2 sets sticky pads + 2-way leads Save \$75 *	FREE		
ENAR Self-Treat Extension Handle	\$145		
ENAR Through-Hair Brush \$120 <input type="checkbox"/> Comb-Brush \$145 <input type="checkbox"/>	\$		
ENAR Probe fine \$100 <input type="checkbox"/> Y large \$145 <input type="checkbox"/> Y small \$100 <input type="checkbox"/>	\$		
ENAR Bell Mushroom Massaging Electrodes	\$175		
ENAR Remote Face Probe <input type="checkbox"/> Acu-Point Probe <input type="checkbox"/>	\$120		
ENAR ET Sturdy 2-Way Leads \$25 <input type="checkbox"/> 4-Way Leads \$70 <input type="checkbox"/>	\$		
ENAR Sticky Pads 8sml <input type="checkbox"/> 4lge <input type="checkbox"/> 2ex-lge <input type="checkbox"/> 2B'fly <input type="checkbox"/>	\$ 25		
Silicon Pads 2 of 75, or 95mm Velcro Band 40/60/80/100 mm	\$ 25		
Conductive Glove \$55 <input type="checkbox"/> Sock \$55 <input type="checkbox"/> Knee \$75 <input type="checkbox"/>	\$		
ENAR Multi-Headed Probes Set of 3 pairs	\$220		
ENAR Therapist Kit incl ENAR, 2 Day Training & Access's Kit.	\$2460		
ENAR Live Training Therapy Basic Certificate - Day 1	\$330		
ENAR Live Training Therapist Advanced Certificate - Day 2	\$330		
ENAR Live Training Integral Beauty Cosmetology - Day 3	\$330		
INFRABED Posture~Pillow Wedge <input type="checkbox"/> EndoBelt/FlexiPad <input type="checkbox"/>	\$440		
INFRABED MiniPad 45cmsq\$440 <input type="checkbox"/> ChairPod 92cm\$880 <input type="checkbox"/>			
INFRABED Medium 150cm \$1360 <input type="checkbox"/> Single 194x80 \$1660 <input type="checkbox"/>			
INFRABED Medium & Single includes BioBlanket Save \$330 <input type="checkbox"/>	FREE		
Energy - BioBlanket Single (120cm x 200cm) 3kg	\$330		
Energy - BioBlanket Double (240cm x 200cm) 6kg	\$660		
Energy - BioBlanket Zip-up Covers (Single \$88 Double \$144)	\$		
+ Freight Costs Devices \$20, Blankets \$25, InfraBeds vary (call us)			
TOTAL FREIGHT International charges vary - (Please request price)			

CREDIT CARD PAYMENTS

Visa Mastercard American Express

Card No. _____ Expiry ___/___ CCV _____

Name on Card _____ Signature _____

BANK TRANSFER to: ANZ Bank, Thirroul Branch - **BSB: 012482 - Account: 198019451**
(SWIFT Code: ANZ BAU 3M) - Identify with your surname as reference, then advise ET.
(Please send COPY of the deposit transaction to info@enlightenedtherapies.com)

ENAR Self-Treat Extension Handle

\$145



The ENAR **Extension Handle** is ideal for self-treating the back and other hard to reach areas. Simply open the arms into their lock out position, clip the ENAR into place and start treating around your torso and over the shoulders.

ENAR Through-Hair Brush (or Comb)

\$120



Through-Hair Brush \$120 (or Comb \$145) ensures excellent skin contact when working directly on the scalp or other areas where there may be abundant hair. Good for cranial neuro stimulation, for stroke, and perfect to treat animals.

ENAR Fine-Point Wand Probe

\$100



This very **fine point beauty electrode** is perfect for treating the face, fine lines, spots and wrinkles. The wand's unique design means it is only for treating yourself, because one end of the circuit is through the hand holding it.

ENAR Bell Mushroom Electrodes

\$175



ENAR **bell mushroom probes** are handheld extension electrodes for massage style remote access. Treat areas large, small or long and patients themselves can deliver treatment to areas that are sensitive or require micro massage.

ENAR Y-Point Probes

\$100



The **Small Y-probe** is an invaluable accessory for treating around toes, fingers, face and neck. The **Large Y-Probe** (\$145) electrode balls are ideal for large and long muscles and for sweeping the spine and the three pathways.

ENAR Acu-Point Probe


\$120



The **Acu-Point Probe** is an extension electrode for treating key-points. ENAR functional for Dose mode, and for personal dental. Acupuncturists use ENAR with the Acu-Point Probe to precisely treat meridians and fine key-points.

ENAR 2-Way Leads


\$25



ENAR **2-Way Leads** are useful for applying ongoing treatment to areas like the back, knees or shoulders, where you can 'set the device and forget it'. With ENAR in a pocket, you can enjoy a treatment whilst moving. New L-Shaped plug.

ENAR 4-Way Leads

\$70



The **4-Way Leads** are able to apply treatment to an even more dispersed area and are useful in energetic techniques where the goal is to lift overall energy. You can enjoy a lengthy treatment while going about other activities.

ENAR Adhesive Pads

\$25



Adhesive Electrode Sticky Pads are reusable. Supplied in a pack of 8 Small pads, 4 Large pads, 2 X-Large or 2 Butterfly pairs. Useful in applying to any areas requiring extended ENAR therapy. Stickiness can be revived with a drop of water.

Long Handle Multi-Headed Probes Set

\$220



This set of long handled **Multi-headed Probes** can be interchanged from convex, to medium or fine headed round probes, simply by unscrewing. Get into narrow areas, or key-points. eg facials.